

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

In re Application of:

Antonius G. P. Oomens et al.

Confirmation No.: 9544

Application Serial No. 10/575,279

Group Art Unit: 1648

Filing Date: April 11, 2006

Examiner: Benjamin P. Blumel

For: RECOMBINANT VIRUSES
WITH HETEROLOGOUS
ENVELOPE PROTEINS

DECLARATION OF ANTONIUS G.P. OOMENS, Ph.D., UNDER 37 C.F.R. § 1.131

Commissioner for Patents
P.O. Box 1450
Alexandria, VA 22313-1450

Sir:

I, Antonius G.P. OOMENS, Ph.D., hereby declare:

1. I am a co-inventor of the above-identified application.
2. I am presenting this Declaration to demonstrate that the claimed subject matter was invented in the U.S. before April, 2003.
3. Attached as Exhibit A is a section of laboratory notebook pages relating to an experiment entitled "Prepare simultaneously 381 and 392 stock to see if GP64 stabilizes RS." On this page, "stock" means a virus suspension; RS = human respiratory syncytial virus (HRSV); 381 = HRSV with glycoproteins deleted and substituted with GP64; 392 = wild-type HRSV with GFP marker gene inserted in place of a non-essential HRSV gene (i.e., control virus).
4. The fourth page of the Exhibit is dated January 3, 2003 ("1/3/03 Timepoint week 8"); the seventh page of the Exhibit is dated, for example, November 20, 2002 ("11/20/02").
5. This experiment was finished in early January, 2003. There is a graph clearly showing that 381, the recombinant virus with substituted GP64, remains stable over an 8-week time period whereas 392, the wild-type control, does not. The straight lines are the main experiment; the dotted lines are a duplication in which the samples were freeze-thawed prior to titrations. The results of the straight and dotted lines were the same, i.e., GP64-containing virus is much more stable. These

notebook pages demonstrate that the subject matter of the above-identified application had been conceived and reduced to practice before April, 2003.

6. I hereby declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true; and further that these statements are made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 of Title 18 of the U.S. Code, and that such willful false statements may jeopardize the validity of the application or any patent issuing thereon.

Date:

1/15/09

signed:




Antonius G.P. Oomens

Exhibit A

11/4 Prepare simultaneously 381 and 392 stock to examine if GP64 stabilizes Ki.

- plated Vero @ 0.8×10^6 / 100 mm dish confluency 90%
~25%
H2 @ 1×10^6 / 100 mm dish ~35%
- 11/5: infect at 16.⁰⁰:
For 392 used 8/20/01 G4 mpt; titer 7.1×10^6 (tubes 1, 2)
For 381 used 8/24/01 G4 mpt; never titered or used before (tubes 1, 2)

For each virus, I infected 3 100mm dishes of H2 and Vero with 10 μ l, 50 μ l, or 250 μ l virus
Virus was added directly to cell medium, and kept on.
Plates were put at 33°C.



- Some percentage infected cells on 11/8, 10.⁰⁰: 60% infected
- Note the decrease of # open cells from 50 \rightarrow 250 μ l for 381
- H2 stock from 381 probably won't work
- For Vero: 381 10 or 50 μ l and 392 10 μ l are most compatible
→ remove 6 out of 10 ml and replace 5 ml fresh DMEM + 5 ml MEM + 5 ml volume now 3 ml, keep at 33°C, (P1, confluency of both cell types ~ 95% at this time)
- Some 11/8 to 11/9
H2 stock: 100% no more virus (100% dead), 100% no more virus (100% dead)

plate 10 and 50 μ l (10) have groups of uninfected healthy He's. Knows like the best amount of virus may be between 50 and 250 μ l.

He/381: Most cells still in good shape, but getting overgrown. Percentage infectivity increasing but too low to get a good titer.

\Rightarrow Since neither 392 nor 381 is perfect, \rightarrow discard He dishes and focus on Vero cells.

Vero/392: 10 μ l appears lowest amount that still leads to near 100% infection. Many syncytia present in the dishes, and medium has turned yellow.

Vero/381: 7 μ l is best dish for some reason as Vero/392 - 10. And because 381 stock may contain many debris or very high particles: pfu ratio judged from the decrease of infectivity in the 250 μ l plate. ^{Medium in 381 dishes was yellow.}

\Rightarrow Choose plates Vero - 10 μ l for both to make stock for exp.

Processing

- Stock's time was on 12/10 at 12:00
- I reduced the sup't from 7.5 \rightarrow 6.2 ml and added 185 μ l 1M hyper \rightarrow bring to 30 mM hyper; color of med. changed to pink. This was done at 13:15 ^{Time \approx 5 days}
- Dishes were further incubated from 13:15 to 17:30, then harvested 2 different ways to ensure that I will get a reasonable stock for both 381 and 392. For each virus:
 1. add another 185 μ l 1M hyper \rightarrow total to \approx 60 mM hyper
 2. scrape cells into medium (\approx 6.5 ml)
 3. PLED 20x ^(and) then divide into 2 15 ml tubes
 4. put one tube at -80°C ^(and) freeze the other as follows:
 - PLED 20x then spin 6 min 750 rpm at RT
 - Xtra cent to resuspend. Resuspend about 100 μ l/tube \rightarrow

~ 8 min in RT H₂O

Σ thaw the frozen tubes and process exactly as the unfrozen one. Store tubes also at 4°C

P5 After aliquotting 100 µl amounts, I changed the order of tubes randomly. I labeled them 381, 381-f/t, 392 & 392-f

Tubes containing freeze-thawed virus are underlined: 381, 392

Titering:

= 4-6 h after putting at 4°C

11/10 For time '0', I took samples from 4°C between 22:00 and 24:00, one at a time, and titered them by TC 1050

1 plated cells at 15:00 earlier in the day; density: 2500 Vfb cells per well in 100 µl

2 set up 7 exp tubes and fill each with 180 µl (one tube will not get virus) DMEM3 + 50 mM hyper

3 add 900 µl DMEM3/hyper to the tube containing 100 µl virus; PUD extensively, then xfer 20 µl to next tube (10^{-2})

4 vortex at setting 5 and repeat the procedure till 10^{-7}

5 Add virus dilutions to 96 well plates:

- vortex again

- electr. pipettor: speed 8, 15 x 10 µl

- use 1st row as uninfected row

Samples titered:

2 x 381

2 x 381-f/t

3 x 392

3 x 392-f/t

11/11 For time '3h', I took samples from 4°C from 13:00 - 16:00 (3 days)

- plate 25 x 10³ Vfb per well in 10 µl in 96 well plate using

- Repeat same as that for timepoint '0' except

• 100 µl virus added to remove contamination from lid

• density: 10⁴ Vfb cells to 8 wells; 100 µl virus to 4 wells; 100 µl

11/18 Timepoint 1 week: done early morning
- protocol exactly as that of time point 3rd.


11/16 Timepoint 2nd week: done early morning. This is actually
a day late, however this should be okay.

12/10 Timepoint week '4'
All (301, 301-1/4, 390, 390-1/11) were done in triplicate this time.

1/3/13 Timepoint week '8'
Same as 12/10.

weeks		392	%	392-f/t	%	381	%	381-f/t	%
0	A	1.20 E7	81.47	3.16 E7	141.51	2.74 E8	100	3.16 E8	61.00
	B	1.61 E7	109.3	2.15 E7	96.28	2.74 E8	100	7.20 E8	139.00
	C	1.61 E7	109.3	1.39 E7	62.25	ND		ND	
average	abc	1.473 E7	100	2.233 E7	100	2.74 E8	100	5.18 E8	100.00
0.5	A	2.42 E7	164.29	1.00 E7	44.78	2.02 E8	73.72	2.42 E8	46.72
	B	1.87 E7	126.95	1.20 E7	53.74	3.16 E8	115.3	4.14 E8	79.92
	C	4.64 E6	31.5	1.00 E7	44.78	ND		ND	
1	A	4.64 E6	31.5	3.57 E6	15.99	3.98 E8	145.26	3.16 E8	61.00
	B	5.35 E6	36.32	3.98 E6	17.82	1.78 E8	64.96	2.51 E8	48.56
	C	5.27 E6	35.78	3.16 E6	14.15	ND		ND	
2	A	3.16 E6	21.45	7.20 E5	3.22	2.51 E8	91.61	3.57 E8	68.92
	B	2.74 E6	18.6	1.39 E6	6.22	3.51 E8	128.1	3.65 E8	70.46
	C	1.78 E6	12.08	2.02 E6	9.05	ND		ND	
4	A	1.39 E5	0.94	1.87 E5	0.84	3.98 E8	145.26	3.51 E8	67.76
	B	1.20 E5	0.81	5.62 E4	0.25	2.74 E8	100	3.16 E8	61.00
	C	5.27 E4	0.36	7.20 E4	0.32	3.16 E8	115.33	3.98 E8	76.84
8	A	3.98 E3	0.03	4.69 E3	0.02	2.02 E8	73.72	4.95 E8	95.56
	B	3.16 E3	0.02	1.61 E3	0.01	3.16 E8	115.33	3.16 E8	61.00
	C	2.80 E3	0.02	3.16 E3	0.01	2.42 E8	88.32	2.42 E8	46.72
12	A								
	B								
	C								
16	A								
	B								
	C								

RSASH vs RSASH,G,F/G^{bac} stability



1	2	3	4	5	6	7	8	9	10	11	12
0 A											
1 B											
2 C											
3 D											
4 E											
5 F	+	+	+	+	+	+	+	+	+	+	+
6 G	+	+	+	+	+	+	+	+	+	+	+
7 H											

3A

3B

1	2	3	4	5	6	7	8	9	10	11	12
0 A											
1 B											
2 C											
3 D											
4 E											
5 F	+	+	+	+	+	+	+	+	+	+	+
6 G	+	+	+	+	+	+	+	+	+	+	+
7 H	+	+	+	+	+	+	+	+	+	+	+

$$\frac{321}{32} = \frac{11/10/02}{11/14/02} - \frac{11/16/02}{11/20/02}$$

321

1	2	3	4	5	6	7	8	9	10	11	12
0 A											
1 B											
2 C											
3 D											
4 E											
5 F											
6 G	+	+	+	+	+	+	+	+	+	+	+
7 H											

3A

3B

1	2	3	4	5	6	7	8	9	10	11	12
0 A											
1 B											
2 C											
3 D											
4 E											
5 F											
6 G	+	+	+	+	+	+	+	+	+	+	+
7 H											

only 3 cells